

## Purpose

The purpose of this lab was to give us a brief introduction as in how materials are moved in and out of cells by mechanisms of passive and active transport. We investigated the basic properties of passive transport like diffusion, osmosis, and differential permeability. Understanding the concept of filtration and the effects of tonicity on cells was also one of the main goals.

## Procedures

### 2-B: Measurement of diffusion through a liquid

1. Working in groups, fill three Petri dishes with 40 ml. of 25°C water.
2. Drop one crystal of potassium permanganate into each dish. Be sure to use the same amount of potassium permanganate for each dish. Record the time.
3. Measure, in millimeters, and record the largest diameter of the colored spot after 5 minutes.
4. Repeat steps 1-3 for water at 5°C and at 45°C.
5. Construct a graph of ranges and means for each temperature.
6. Based on your knowledge of diffusion, what is an explanation for these results?

### 2-C: Measurement of diffusion through agar

1. Petri dishes have been filled with agar. Two holes have been made in the agar. Into one hole, place two drops of methylene blue. Into the other hole, place two drops of potassium permanganate. Record the time and immediate diameter of each spot. This will be your time zero measurement.
2. Measure the diameter of each spot in millimeters once every minute for fifteen minutes. Calculate the averages from the data collected by all groups doing this exercise. Summarize these data.
3. Construct a graph of average diffusion diameter versus time for both chemicals.
4. Determine the diffusion rate of each chemical. Which has the fastest diffusion rate, methylene blue or potassium permanganate? Record these results.
5. Look up the molecular formula and structure of methylene blue and potassium permanganate in a Merck Index. Make note of this information.
6. Interpret your result with respect to the information obtained from the Merck Index.

### 2-D: Demonstration of filtration

1. Fold three filter papers into cones and insert them into three separate glass funnels. Wet

the papers to make them stick to the glass.

2. Prepare three 100-milliliter solutions of charcoal and water. Make one thick, one medium thickness, and one thin. Record the mass of the charcoal used in each preparation.

NOTE: if your "thin" solution continually runs through the filter, making it impossible to count drops, it is too thin; you will need to make all your solutions proportionally thicker.

3. Pour 50 ml of each solution, one at a time, into a funnel.

4. Immediately count the number of drops produced per minute.

NOTE: it may be easier to count the drops for 15 seconds then multiply by four to obtain drops per minute.

5. Count the number of drops per minute when the funnel is half-filled.

6. Count the number of drops per minute when the funnel is nearly empty.

7. Did the charcoal pass into the filtrate? Which solution had the fastest rate of filtration? What is the driving force behind filtration? What other factors influence the rate of filtration? Do your results illustrate these influencing factors?

8. Repeat these procedures with the remaining 50 ml. of solution.

## 2-F: Measurement of osmosis

1. Attach dialysis bags filled as much as possible with sucrose solutions securely to the bottom of two open, thin glass tubes. One bag should be filled with a 25% sucrose solution and the other should be filled with a 50% sucrose solution. Make sure ends of the tubes are immersed in the solutions.

NOTE: reliable results depend on your ability to tightly seal the dialysis bags.

2. Insert both bags into separate beakers of distilled water making sure the dialysis bags are fully submersed but not touching the bottom of the beakers, and suspend each by gently applying a ring stand clamp to the glass tubes. Check for solution leaking out of the bags.

3. Allow five minutes for the systems to equilibrate. Then, mark the fluid levels of each glass tube with a felt pen. Record the time.

4. Weight the glass tubes.

5. If the fluid level rises to the top of the glass tube sooner than 50 minutes, record the time it took to get there, measure the length in millimeters from the equilibration line to

the top of glass tube. Divide that length by the number of minutes to get your rate in mm/min.

6. Determine the rate of osmosis for each system. Which system had the fastest osmotic rate, the 25% or 50% sucrose solution? Explain these results.

## 2-G: Measurement of differential permeability of sugar and starch

### Procedure

1. Fill a dialysis bag with a 1% starch – 10% glucose solution. Reliable results depend on

your ability to tightly seal the dialysis bag.

2. Tie the bag to a glass rod and suspend it in a beaker of distilled water.

NOTE: Test the water from the bottom of the beaker to ensure that it is free of starch and/or sugar.

3. After 15 minutes has passed check the water again for starch and sugar in the following way:

Test for starch:

a. Add 10 drops of Lugol's solution to 5 ml of water obtained from the beaker.

Reddish color = No starch

Navy blue color = Starch present

Test for sugar:

a. Add 3 ml of Benedict's solution to 5 ml of water obtained from the beaker.

Simmer the solution at a low boil for 5 minutes.

Blue color = No sugar

Color change = Sugar present (green = little sugar; yellow = moderate sugar;

orange = more sugar; red = lots of sugar) 4. Test the water in the beaker again at 30, 45 and 60 minutes.

5. Record these results. Explain the significance of these findings in relation to the permeability of the dialysis bag.

1. One milliliter of each of the following solutions will be in three separate test tubes.

a. Distilled water (hypotonic)


b. Physiological saline – 0.85% NaCl (isotonic) c. Salt water – 2.0% NaCl (hypertonic)

## 2-H: The effects of tonicity on red blood cells - Demonstration


2. A small drop of blood will be added to each tube and the contents thoroughly mixed.
3. A wet mount slide will be made of each solution.
4. Examine each slide under the high-dry lens of a compound microscope.
5. Observe the following:
  - a. Hemolysis of cells in the hypotonic solution. (Note the transparent solution.)
  - b. Maintenance of cell size in the isotonic solution.
  - c. Crenation of cells in the hypertonic solution.
6. Make a drawing of each observation and provide an explanation for each.

## Results

### 2-B: Measurement of diffusion through a liquid



### 2-C: Measurement of diffusion through agar




## 2-G: Measurement of differential permeability of sugar and starch

| Time(min) | Lugol's | Benedict's | Interpretation       |
|-----------|---------|------------|----------------------|
| 15        | Reddish | Blue       | N/A present          |
| 30        | Reddish | Blue       | N/A present          |
| 45        | Reddish | Green/Blue | Slight sugar present |
| 60        | Reddish | Green      | Slight sugar present |

Since there was a color change to green when Benedict's solution was added to the 5ml of water, it indicated that there was a little bit of sugar present inside the dialysis bag. These results were due to my team not being able to successfully tie the bag tight enough.

## 2-H: The effects of tonicity on red blood cells – Demonstration



A. Cells in a hypertonic solution will shrivel and die because through plasmolysis.

B. Cells in distilled water, a hypotonic solution, will swell and expand until they burst because of cytolysis.

C. Cells in an isotonic solution will be happy because the salt and water is balanced, nothing happens.

## **Discussion**

I think this was my least favorite lab because there was just too many different procedures to go complete and even though they all had a similar purpose, I did not

think all of them were necessary. Some took a long time, and I didn't really get the point of doing them. The only interesting thing about this lab was learning how blood cells can't survive with in different solutions because I felt like that actually had to do with real life.

## **Conclusion**

- Understand the mechanism of Brownian motion.
- Understand the difference between passive and active transport.
- Be able to define diffusion, osmosis, active transport, dialysis, and filtration.
- Know the result of dropping red blood cells in hypertonic, isotonic, and hypotonic solutions.
- Understand the significance of all of these experiments in terms of passive transport processes and molecular activity.